


0.01
ct and above




Initially all the goods of size 0.01 ct and above are tested on **Diamond Sure or AMS**



All the stones passed in above machines are mixed in the stock while the referred are tested further on **Diamond View**



All the stones passed on Diamond View are mixed in stock while the referred stones are further tested on **Enspectre R532 (Raman)**



All the stones passed on Enspectre R532(Raman) are mixed in stock while the failed stones (synthetic) are sent to an **External Lab** for double confirmation

0.01
ct and below



The stones are initially screened on **Q-CHK**, the passed stones are mixed in stock while the synthetic stones as indicated by the machine are further tested on Enspectre-532 (Raman)



The stones passed on Enspectre R532(Raman) are mixed in stock while the failed stones (synthetic) are sent to an **External Lab** for double confirmation